


Processing Method:	Derivati PhNCO_2	Project Name:	Interesterifikacija\Deriv PhNCO
Processing Method ID:	1545	System:	Massspek
Channel:	234.0nm	Calibration ID:	1546
Proc. Chnl. Descr.:	****		
Date Calibrated:	1/23/2019 5:06:48 PM EET		

A 1.821973e+005
 B 1.376691e+009
 C 0.000000e+000
 D 0.000000e+000
 R^2 0.999721

Calibration Plot



Peak: Glicerins

	Name	Level	X Value	Response	Calc. Value	% Deviation	Manual	Ignore
1	Glicerins	2-100	0.000239	384313.307020	0.000147	-38.572	Yes	No
2	Glicerins	2-10	0.002390	3627957.277926	0.002503	4.725	Yes	No
3	Glicerins	2-0.5	0.011950	16605098.440474	0.011929	-0.174	Yes	No
4	Glicerins	2-1	0.023901	22340389.627700			Yes	Yes